


GreenGel™ with GelGreen™

GG3600 1% Agarose TAE w/ GelGreen, Landscape Gel 12+1 well



Run @ 100 V for 45 min.

Lane# 1 & 12: SmartMark 1 Kb DNA ladder, 40 ng per band
Lane# 2 & 11: Lambda/EcoRI, top band is 200 ng
Lane# 3, 6 & 9: Plasmid DNA
Lane# 4, 7 & 10: 3 Kb PCR fragment, 300 ng ea.
Lane# 5: 5 Kb PCR fragment, 100 ng ea.
Lane# 8: 2112 bp PCR Fragment



No HazMat

Quick and easy disposal in your local trash can.
You don't always have to look like an astronaut to throw stuff away...



No UV necessary

Use blue light instead and eliminate goggles and masks.
You can put away the beach towel and tanning lotion as well!



No Photo-nicking

Better DNA for downstream application.
Nobody likes their photos damaged or nicked.



No Half & Half

Uniform brightness rather than half dark, half light as with EtBr gels.
If you like your coffee black you'll love these gels!

GreenGel with GelGreen

Catalog #	Gel Size (cm)	Description	List Price 1 bx	20% Discount 2-5 bx	30% Discount 6+ bx
GG5100	Portrait 6.0(W) x 9.5 (L)	1% Agarose TAE w/GelGreen, Portrait Gel, 5+1 well, 1.5mm comb, 10/bx (5 wells - 45µl and a half size reference well for standard)	\$97.95	\$78.36	\$68.57
GG3870	Portrait 6.0(W) x 9.5 (L)	1% Agarose TAE w/GelGreen, Portrait Gel, 10 well, 1.0mm comb, 10/bx (10 wells - 15µl per well)	\$97.95	\$78.36	\$68.57
GG3600	Landscape 11.5(W) x 5.6 (L)	1% Agarose TAE w/GelGreen, Landscape Gel, 12+1 well, 1.0mm comb, 10/bx (12 wells - 30µl and a smaller size reference well for standard)	\$97.95	\$78.36	\$68.57
GG3750	Landscape 11.5(W) x 5.6 (L)	1% Agarose TAE w/GelGreen, Landscape Gel, 24+1 well, 1.0mm comb, 10/bx (24 wells - 11µl and a smaller size reference well for standard)	\$97.95	\$78.36	\$68.57
GG3690	Medium 11.5(W) x 8.9 (L)	1% Agarose TAE w/GelGreen, Medium Gel, 5+2 well, 1.5mm comb, 5/bx (2 wells - 112µl, 3 wells - 72µl and 2 smaller size wells for standard)	\$97.95	\$78.36	\$68.57
GG9999	Mix & Match	Sample Package of 1% Agarose TAE w/ GelGreen: mix & match of any 3 packs of Portrait, Landscape or Medium Gel; 3 packs per box	\$49.95	N/A	N/A

Promotion for standard gels listed above only. Promotion ends December 31, 2014