





7738 Arjons Drive
 San Diego, CA 92126
 858 684-3190
 inquiry@embitec.com

Presorted Standard
 U.S. Postage
 PAID
 San Diego, CA
 Permit #2686

The Light Source is the Key
 Make Your Samples Easy to See

ViewOne LabLite Set

ViewOne™ LabLite Set




*Includes ViewOne LabLite and all accessories

Catalog No.
 VL-1200
 \$199



ViewOne LabLite



ViewOne to View All!

- It may be tiny, but it's versatile - use it to view your Coomassie Blue gels, help you pipette your plates or pick colonies!
- And that's all in a sleek, water-resistant enclosure, thinner than the standard CD case!

When things are small and hard to see,
a ViewOne is all you need - take it from me!

-Quote from old guy in lab

Catalog No.
VL-1000
\$145

But wait, there are accessories...

Prop It Up with Stands!

- Work on the ViewOne LabLite at a 15° or 25° angle.
- Propped high enough for the comfortable viewing angle, without compromising the comfortable pipetting angle.



96-well Template!

- Go ahead, mark and label what you need without intrusively marking up your plate.
- Now there's no confusion about which well you're picking.

